

ANTERIOR TEETH


SR **PHONARES**[®]
NHC

SOFT










YOUTHFUL

UNIVERSAL

MATURE



BOLD

	YOUTHFUL	UNIVERSAL	MATURE
Small	<p>B61</p>  <p>42.2</p> <p>10.7</p> <p>8.2</p>	<p>B71</p>  <p>43.2</p> <p>10.3</p> <p>7.9</p>	<p>B81</p>  <p>42.0</p> <p>9.7</p> <p>7.6</p>
Medium	<p>B62</p>  <p>45.8</p> <p>11.8</p> <p>9.0</p>	<p>B72</p>  <p>45.0</p> <p>10.8</p> <p>8.3</p>	<p>B82</p>  <p>45.6</p> <p>10.7</p> <p>8.3</p>
Large	<p>B63</p>  <p>50.2</p> <p>13.2</p> <p>9.8</p>	<p>B73</p>  <p>49.6</p> <p>11.9</p> <p>9.0</p>	<p>B83</p>  <p>49.0</p> <p>11.1</p> <p>9.2</p>




POSTERIOR TEETH

SR PHONARES[®]
Typ NHC



SR PHONARES[®]
Lingual NHC



		SR PHONARES NHC		SR PHONARES <i>Typ</i> NHC	SR PHONARES <i>Lingual</i> NHC
		upper	lower		
Small	soft	S61	L50	NU3 / NL3	LU3 / LL3
		S71	L51	NU3 / NL3, NU5 / NL5	LU3 / LL3, LU5 / LL5
		S81	L50, L51	NU3 / NL3	LU3 / LL3, LU5 / LL5
	bold	B61	L50, L51	NU3 / NL3	LU3 / LL3
		B71	L50, L51	NU3 / NL3, NU5 / NL5	LU3 / LL3
		B81	L50, L51	NU3 / NL3	LU3 / LL3
Medium	soft	S62	L51	NU3 / NL3, NU5 / NL5	LU5 / LL5
		S72	L51, L53	NU5 / NL5	LU5 / LL5, LU6 / LL6
		S82	L51, L53	NU5 / NL5	LU5 / LL5, LU6 / LL6
	bold	B62	L51, L53	NU3 / NL3, NU5 / NL5	LU5 / LL5, LU6 / LL6
		B72	L51	NU3 / NL3, NU5 / NL5	LU5 / LL5, LU6 / LL6
		B82	L51, L53	NU3 / NL3, NU5 / NL5	LU3 / LL3, LU5 / LL5
Large	soft	S63	L52, L53	NU5 / NL5	LU6 / LL6
		S73	L52	NU5 / NL5, NU6 / NL6	LU6 / LL6
		S83	L52, L53	NU5 / NL5, NU6 / NL6	LU6 / LL6
	bold	B63	L52, L53	NU6 / NL6	LU6 / LL6
		B73	L52, L53	NU6 / NL6	LU6 / LL6
		B83	L52	NU5 / NL5, NU6 / NL6	LU6 / LL6

This combination table is recommended as a guideline. In case of particular anatomical conditions, deviations are possible.